

INTERACRYL PLAST

IA/2463

- (SL) Navodila za uporabo / (EN) Directions for use
(DE) Gebrauchsanweisung / (HR) Uputa za korištenje
(RU) Инструкции по использованию / (MK) Препорака за употреба
(CS) Návod na zpracování / (ES) Instrucciones de uso
(KK) Пайдалану бойынша нұсқаулық / (GR) Διαδικασία επεξεργασίας
(RO) Instrucțiuni de utilizare / (ALB) Pergatitja

CS

INTERDENT s.r.o.
Foerstrova 12, Strašnice
CZ-10000 Praga
T: +420/274 783 114
F: +420/274 820 130
E: interdent@interdent.cz

HR

INTERDENT d.o.o
Vinogradski odvojak 2d
HR-10431 Sveta Nedelja
T: +385/1 3873 644
F: +385/1 38736 17
E: interdent@interdent.hr

SK

INTERDENT SK S.R.O.
Za dráhou 21
SK-902 01 Pezinok
T: 0903 418 001
E: interdent@interdent.sk

SR

INTERDENT d.o.o.
Egipatska 22, lok 3
RS-11070 Novi Beograd
T/F: +381/11 217 53 74
www.interdent-bg.com



 **ALPHABOND** dental

Ph: (02) 9417 6660 www.alphabond.com.au

Verzija: 08/2023

Datum: 18.08.2023

Made in Slovenia



Potrebno upoštevati navodila za uporabo / Consult instructions for use / Bitte beachten Sie die Bedienungsanleitung/Consultar las instrucciones de uso/Potrebno poštivati uputstva za upotrebu/Водете се по упатството за употреба / Необходимо учитывать инструкции по применению / Viz návod k použití / Қолдану бойынша нұсқаулықты ескеру қажет / Consultati instructiunile de folosire / Konsultohuni me udhëzimet për përdorim / Συμβουλευτείτε τις οδηγίες χρήσης



Uporabno do / Use by / Gültig bis / Usar antes / Upotrijebiti do / Рок на употреба до / Использовать до / Použítí / Дейін қолданыңыз / Folosit de catre / Përdorimi sipas / Χρήση με



Kontrolna številka / Batch code / Seriennummer / Numero de serie / Serijski broj / Серијски број / Серийный номер / Seriové číslo / Сериялық нөмір / Codul lotului / Kodi i grupit (grumbullit) / Κωδικός παρτίδας



Kataloška številka / Catalogue number / Katalognummer / Numero de catalogo / Kataloški broj / Каталогски број / Каталогный номер / Katalogové číslo / Каталогтық нөмір / Numarul de catalog/Numri i katalogut / Αριθμός καταλόγου



Medicinski pripomoček / Medical device / Medizinprodukt / Medicinski proizvod / / Медицинское изделие / Медицински уред / Zdravotnický prostředek / Productio sanitario / Ιατρική συσκευή / Медициналық құрылғы / Dispozitiv medical / Pajisje mjekësore



Edinstveni identifikator pripomočka / Unique device identifier / Einmalige Produktkennung / Jedinstvena identifikacija proizvoda / / Уникальный идентификатор изделия / Единствен идентификатор на уред / Jedinečným identifikátorem prostředku / Identificador único del producto / Μοναδικό αναγνωριστικό συσκευής / Бірегей өнім идентификаторы / Identificator unic de dispozitiv / Identifikues unik i pajisjes



Proizvajalec / Manufacturer / Hersteller / Fabricante / Proizvođač / Производител / Производител / Výrobce / Өндіруші / Producator / Prodhuiesi / Κατασκευαστής



Datum proizvodnje / Date of manufacture / Herstellungsdatum / Fecha de fabricacion / Datum proizvodnje / Датум на производство / Дата изготовления / Datum výroby / Өндірілген күні / Data fabricatiei / Data e prodhimit / Ημερομηνία κατασκευής



Hraniti zaščiteno pred sončno svetlobo / Keep away from sunlight / Geschützt gegen Sonnenlicht halten / Čuvati odvojeno od sunčeve svjetlosti / Хранение -вдали от солнечного света / Чувајте го подалеку од сончева светлина / Uchovávejte mimo dosah slunečního záření / Mantener alejado de la luz solar / Күн жарығынан ұзақ сақтаңыз / A se tine departe de lumina soarelui / Mbani larg nga rrezet e diellit / Κρατήστε το μακριά από ηλιακή ακτινοβολία



Lomljivo, ravhati previdno / Fragile, handle with care / Zerbrechlich, Vorsichtig behandeln / Lomljivo, pažljivo rukovati / Křehké, zacházejte opatrně / Frágil, manipular con cuidado / Хрупкое, обращаться с осторожностью / Кревки, внимателно се справувајте / Сыңғыш, байқап ұстаныз / Εύθραυστο, χειρισμός με προσοχή / Fragil, utilizati cu grija / E brishtë, trajto me kujdes.

INTERACRYL PLAST

(Acrylic material for custom impression trays)

Interacryl plast is a cadmium-free, self-curing material based on methyl methacrylate. It is supplied as a polymer and monomer and is intended for fabricating custom impression trays.

Indications:

- Fabricating custom impression trays.

Contraindications:

- not for removable dentures
- not for orthodontic

Directions for use:**Mixing ratio at room temperature:**

Mixing ratio: 33,6 g powder / 8 ml liquid (4,2 g powder / 1 ml liquid)

First pour the liquid into a silicone mixing vessel, then add the required amount of powder. Mix the two components, using a spatula, until the dough is no longer sticky and separates cleanly from the sides of the mixing vessel. This should take about 90 seconds at room temperature (22°C). Then the dough can be applied directly to the plaster model - which should have been coated with Isoacryl in advance. During the plastic phase, the material can be processed for approximately 6 minutes. At room temperature, the tray material will cure 12 minutes after mixing. To ensure that the tray retains its shape optimally, it should be post-cured by immersing it in water at a temperature of 60°C for a period of 4 minutes.

Caution:

Do not mix polymer powder with another supplier monomer liquid and vice versa.

Danger: Highly flammable liquid and vapour. / May cause respiratory irritation. / Causes skin irritation. / May cause an allergic skin reaction. / Causes serious eye irritation. / Keep away from heat/sparks/open flame/hot surfaces. - No smoking. / Keep container tightly closed. / Avoid breathing vapours. / Keep away from heat/sparks/open flames/hot surfaces. - No smoking / Keep container tightly closed. Avoid breathing vapours. / IF ON SKIN: Wash with plenty of soap and water. / Wear protective gloves/protective clothing/eye protection/face protection. / IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing./ If eye irritation persists: Get medical advice/attention. / Dispose of contents/containers in accordance with local regulation.

Personal protection:

We recommend wearing half facemask with filter A1.

Storage:

- at room temperature
- in the place with good ventilation
- keep away from sunlight
- keep away from children

Warranty:

Any technical consulting regarding our products, whether verbal, in writing or during demonstrations, is carried out to the best of our knowledge and is to be considered nonbinding. We have no control over use, application and processing, for which the user is responsible.